



TM4SF5 Antibody, HRP conjugated

Product Code	CSB-PA023620LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14894
Immunogen	Recombinant Human Transmembrane 4 L6 family member 5 protein (112-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane 4 L6 family member 5 (Tetraspan transmembrane protein L6H), TM4SF5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TM4SF5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.