



TMEM140 Antibody, Biotin conjugated

Product Code	CSB-PA023710LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NV12
Immunogen	Recombinant Human Transmembrane protein 140 protein (33-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 140, TMEM140
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM140
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.