



TMEM184C Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA023771LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9NVA4 |
| Immunogen | Recombinant Human Transmembrane protein 184C protein (107-178AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transmembrane protein 184C (Transmembrane protein 34), TMEM184C, TMEM34 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | TMEM184C |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |