



TMEM196 Antibody, Biotin conjugated

Product Code	CSB-PA023788LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HYL7
Immunogen	Recombinant Human Transmembrane protein 196 protein (53-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 196, TMEM196
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM196
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.