





TMEM31 Antibody

Product Code	CSB-PA023828LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JXX7
Immunogen	Recombinant Human Transmembrane protein 31 protein (1-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 31, TMEM31
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM31


Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA023828LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA023828LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA023828LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.