



TMEM47 Antibody, FITC conjugated

Product Code	CSB-PA023845LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQJ4
Immunogen	Recombinant Human Transmembrane protein 47 protein (42-82AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 47 (Brain cell membrane protein 1) (Transmembrane 4 superfamily member 10), TMEM47, BCMP1 TM4SF10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM47
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.