



TMEM68 Antibody, Biotin conjugated

Product Code	CSB-PA023872LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MH6
Immunogen	Recombinant Human Transmembrane protein 68 protein (74-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 68, TMEM68
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMEM68
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.