




TMEM91 Antibody


| | |
|----------------------------|---|
| Product Code | CSB-PA023896LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6ZNR0 |
| Immunogen | Recombinant Human Transmembrane protein 91 protein (1-97AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB, IHC, IF; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IF:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transmembrane protein 91 (Dispanin subfamily C member 3) (DSPC3), TMEM91 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | TMEM91 |

Image




Western Blot


All lanes: TMEM91 antibody at 16µg/ml + A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 19, 16, 14, 15 kDa
 Observed band size: 19, 34 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA023896LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA023896LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA023896LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.