




TMEM91 Antibody

Product Code	CSB-PA023896LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZNR0
Immunogen	Recombinant Human Transmembrane protein 91 protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 91 (Dispanin subfamily C member 3) (DSPC3), TMEM91
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM91

Image



Western Blot


All lanes: TMEM91 antibody at 16µg/ml + A549 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 19, 16, 14, 15 kDa


Observed band size: 19, 34 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA023896LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA023896LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA023896LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)