



TMSB4Y Antibody

Product Code	CSB-PA023935LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14604
Immunogen	Recombinant Human Thymosin beta-4, Y-chromosomal protein (2-44AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymosin beta-4, Y-chromosomal, TMSB4Y, TB4Y
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMSB4Y