



Tnfsf11 Antibody, HRP conjugated

Product Code	CSB-PA023986LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35235
Immunogen	Recombinant Mouse Tumor necrosis factor ligand superfamily member 11 protein (73-306AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 11 (Osteoclast differentiation factor) (ODF) (Osteoprotegerin ligand) (OPGL) (Receptor activator of nuclear factor kappa-B ligand) (RANKL) (TNF-related activation-induced cytokine) (TRANCE) (CD antigen CD254) [Cleaved into: Tumor necrosis factor ligand superfamily member 11, membrane form; Tumor necrosis factor ligand superfamily member 11, soluble form], Tnfsf11, Opgl Rankl Trance
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Tnfsf11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.