





TPSB2 Antibody

Product Code	CSB-PA024128LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20231
Immunogen	Recombinant Human Tryptase beta-2 protein (31-275AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptase beta-2 (Tryptase-2) (EC 3.4.21.59) (Tryptase II), TPSB2, TPS2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TPSB2

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA024128LA01HU at dilution of 1:100

**Western Blot**

Positive WB detected in: Rat liver tissue, Rat kidney tissue, Mouse liver tissue, Mouse kidney tissue, Mouse stomach tissue

All lanes: TPSB2 antibody at 3 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.