



TPSB2 Antibody, Biotin conjugated

Product Code	CSB-PA024128LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20231
Immunogen	Recombinant Human Tryptase beta-2 protein (31-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptase beta-2 (Tryptase-2) (EC 3.4.21.59) (Tryptase II), TPSB2, TPS2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TPSB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.