




TSLP Antibody


Product Code	CSB-PA025141LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969D9
Immunogen	Recombinant Human Thymic stromal lymphopoietin protein (29-159AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymic stromal lymphopoietin, TSLP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TSLP

Image



Western Blot

Positive WB detected in: U87 whole cell lysate, SH-SY5Y whole cell lysate, Mouse liver tissue
 All lanes: TSLP antibody at 3.4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 19, 8 kDa
 Observed band size: 19 kDa



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA025141LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.