



TSPAN8 Antibody, Biotin conjugated

Product Code	CSB-PA025166LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19075
Immunogen	Recombinant Human Tetraspanin-8 protein (110-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-8 (Tspan-8) (Transmembrane 4 superfamily member 3) (Tumor-associated antigen CO-029), TSPAN8, TM4SF3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TSPAN8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.