



UBE2V1 Antibody, FITC conjugated

Product Code	CSB-PA025482LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13404
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 variant 1 protein (2-147AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 variant 1 (UEV-1) (CROC-1) (TRAF6-regulated IKK activator 1 beta Uev1A), UBE2V1, CROC1 UBE2V UEV1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UBE2V1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.