




UBL3 Antibody

Product Code	CSB-PA025498LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95164
Immunogen	Recombinant Human Ubiquitin-like protein 3 protein (17-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like protein 3 (Membrane-anchored ubiquitin-fold protein) (HsMUB) (MUB) (Protein HCG-1), UBL3, PNSC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UBL3

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA025498LA01HU at dilution of 1:100