



UBL3 Antibody, FITC conjugated

Product Code	CSB-PA025498LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95164
Immunogen	Recombinant Human Ubiquitin-like protein 3 protein (17-114AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like protein 3 (Membrane-anchored ubiquitin-fold protein) (HsMUB) (MUB) (Protein HCG-1), UBL3, PNSC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UBL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.