





UBXN8 Antibody

Product Code	CSB-PA025539LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00124
Immunogen	Recombinant Human UBX domain-containing protein 8 protein (55-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UBX domain-containing protein 8 (Reproduction 8 protein) (Rep-8 protein) (UBX domain-containing protein 6), UBXN8, D8S2298E REP8 UBXD6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBXN8

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA025539LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA025539LA01HU at dilution of 1:100