



UBXN8 Antibody, HRP conjugated

Product Code	CSB-PA025539LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00124
Immunogen	Recombinant Human UBX domain-containing protein 8 protein (55-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UBX domain-containing protein 8 (Reproduction 8 protein) (Rep-8 protein) (UBX domain-containing protein 6), UBXN8, D8S2298E REP8 UBXD6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBXN8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.