



UCP2 Antibody, HRP conjugated

Product Code	CSB-PA025555EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55851
Immunogen	Recombinant Human Mitochondrial uncoupling protein 2 protein (137-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial uncoupling protein 2 (UCP 2) (Solute carrier family 25 member 8) (UCPH), UCP2, SLC25A8
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	UCP2