



UCP2 Antibody, Biotin conjugated

Product Code	CSB-PA025555ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55851
Immunogen	Recombinant Human Mitochondrial uncoupling protein 2 protein (137-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial uncoupling protein 2 (UCP 2) (Solute carrier family 25 member 8) (UCPH), UCP2, SLC25A8
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	UCP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.