



UCP3 Antibody, Biotin conjugated

Product Code	CSB-PA025556LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55916
Immunogen	Recombinant Human Mitochondrial uncoupling protein 3 protein (33-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial uncoupling protein 3 (UCP 3) (Solute carrier family 25 member 9), UCP3, SLC25A9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UCP3