



UROS Antibody, Biotin conjugated

Product Code	CSB-PA025679LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10746
Immunogen	Recombinant Human Uroporphyrinogen-III synthase protein (1-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uroporphyrinogen-III synthase (UROIIIIS) (UROS) (EC 4.2.1.75) (Hydroxymethylbilane hydrolyase [cyclizing]) (Uroporphyrinogen-III cosynthase), UROS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UROS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.