



USP5 Antibody, FITC conjugated

Product Code	CSB-PA025742LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45974
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 5 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 5 (EC 3.4.19.12) (Deubiquitinating enzyme 5) (Isopeptidase T) (Ubiquitin thioesterase 5) (Ubiquitin-specific-processing protease 5), USP5, ISOT
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.