



VAMP1 Antibody, FITC conjugated

Product Code	CSB-PA025780LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23763
Immunogen	Recombinant Human Vesicle-associated membrane protein 1 protein (1-96AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle-associated membrane protein 1 (VAMP-1) (Synaptobrevin-1), VAMP1, SYB1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VAMP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.