



VAMP5 Antibody, HRP conjugated

Product Code	CSB-PA025784LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95183
Immunogen	Recombinant Human Vesicle-associated membrane protein 5 protein (1-72AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle-associated membrane protein 5 (VAMP-5) (Myobrevin), VAMP5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VAMP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.