



WDR62 Antibody, FITC conjugated

Product Code	CSB-PA026050LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43379
Immunogen	Recombinant Human WD repeat-containing protein 62 protein (889-1188AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 62, WDR62, C19orf14
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	WDR62
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.