




XG Antibody

Product Code	CSB-PA026190LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55808
Immunogen	Recombinant Human Glycoprotein Xg protein (22-143AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycoprotein Xg (Protein PBDX), XG, PBDX
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	XG

Image



Western blot

All lanes: XG antibody at 10µg/ml + Mouse heart tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 20, 22 kDa

Observed band size: 20, 37, 60 kDa