



XG Antibody, HRP conjugated

Product Code	CSB-PA026190LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55808
Immunogen	Recombinant Human Glycoprotein Xg protein (22-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycoprotein Xg (Protein PBDX), XG, PBDX
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	XG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.