




ZNF215 Antibody

Product Code	CSB-PA026585LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL58
Immunogen	Recombinant Human Zinc finger protein 215 protein (12-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 215 (BWSCR2-associated zinc finger protein 2) (BAZ-2) (Zinc finger protein with KRAB and SCAN domains 11), ZNF215, BAZ2 ZKSCAN11
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF215

Image



Western blot

All lanes: ZNF215 antibody at 3µg/ml

Lane 1: Hela whole cell lysate


Lane 2: A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 61, 36 kDa

Observed band size: 61 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA026585LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.