



ZNF215 Antibody, FITC conjugated

Product Code	CSB-PA026585LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL58
Immunogen	Recombinant Human Zinc finger protein 215 protein (12-304AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 215 (BWSCR2-associated zinc finger protein 2) (BAZ-2) (Zinc finger protein with KRAB and SCAN domains 11), ZNF215, BAZ2 ZKSCAN11
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF215