



ZNF215 Antibody, Biotin conjugated

Product Code	CSB-PA026585LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL58
Immunogen	Recombinant Human Zinc finger protein 215 protein (12-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 215 (BWSCR2-associated zinc finger protein 2) (BAZ-2) (Zinc finger protein with KRAB and SCAN domains 11), ZNF215, BAZ2 ZKSCAN11
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF215
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.