



ZNF366 Antibody, HRP conjugated

Product Code	CSB-PA026702LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N895
Immunogen	Recombinant Human Zinc finger protein 366 protein (455-650AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 366 (Dendritic cell-specific transcript protein) (DC-SCRIPT), ZNF366
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF366
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.