



# ZNF385C Antibody, HRP conjugated

Product Code	CSB-PA026712LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q66K41
Immunogen	Recombinant Human Zinc finger protein 385C protein (1-174AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 385C, ZNF385C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF385C