



ZNF385C Antibody, FITC conjugated

Product Code	CSB-PA026712LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q66K41
Immunogen	Recombinant Human Zinc finger protein 385C protein (1-174AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 385C, ZNF385C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF385C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.