



ZNF398 Antibody, Biotin conjugated

Product Code	CSB-PA026722LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD17
Immunogen	Recombinant Human Zinc finger protein 398 protein (401-642AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 398 (Zinc finger DNA-binding protein p52/p71), ZNF398, KIAA1339 ZER6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF398
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.