



ZNF410 Antibody, Biotin conjugated

Product Code	CSB-PA026730LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VK4
Immunogen	Recombinant Human Zinc finger protein 410 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 410 (Another partner for ARF 1) (Zinc finger protein APA-1), ZNF410, APA1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF410
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.