



ZNF486 Antibody, Biotin conjugated

Product Code	CSB-PA026783LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96H40
Immunogen	Recombinant Human Zinc finger protein 486 protein (251-463AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 486 (KRAB domain only protein 2), ZNF486, KRBO2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF486
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.