





# ZNF529 Antibody

<b>Product Code</b>	CSB-PA026817LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6P280
<b>Immunogen</b>	Recombinant Human Zinc finger protein 529 protein (393-563AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Zinc finger protein 529, ZNF529, KIAA1615
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	ZNF529


## Image




Western blot  
 All lanes: ZNF529 antibody at 6µg/ml  
 Lane 1: Hela whole cell lysate  
 Lane 2: MCF-7 whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 66 kDa  
 Observed band size: 66 kDa



Immunohistochemistry analysis of human lung cancer using CSB-PA026817LA01HU at dilution of 1:100



Immunohistochemistry analysis of human glioma using CSB-PA026817LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA026817LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.