



ZNF530 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA026818LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6P9A1 |
| Immunogen | Recombinant Human Zinc finger protein 530 protein (301-599AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc finger protein 530, ZNF530, KIAA1508 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | ZNF530 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |