




ZNF554 Antibody

Product Code	CSB-PA026834LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TJ5
Immunogen	Recombinant Human Zinc finger protein 554 protein (239-538AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 554, ZNF554
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF554

Image



Western blot

All lanes: ZNF554 antibody at 5µg/ml

Lane 1: K562 whole cell lysate

Lane 2: U937 whole cell lysate

Lane 3: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 61 kDa

Observed band size: 61 kDa