



ZNF586 Antibody, Biotin conjugated

Product Code	CSB-PA026869LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXT0
Immunogen	Recombinant Human Zinc finger protein 586 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 586, ZNF586
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF586
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.