



ZNF627 Antibody, HRP conjugated

Product Code	CSB-PA026905LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L945
Immunogen	Recombinant Human Zinc finger protein 627 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 627, ZNF627
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF627
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.