



ZNF765 Antibody

Product Code	CSB-PA027021LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L2R6
Immunogen	Recombinant Human Zinc finger protein 765 protein (1-60AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 765, ZNF765
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF765
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.