




ZNF821 Antibody

Product Code	CSB-PA027066LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75541
Immunogen	Recombinant Human Zinc finger protein 821 protein (1-370AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 821, ZNF821
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF821

Image



Western blot


All lanes: ZNF821 antibody at 8µg/ml + MCF-7 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 47, 43 kDa

Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA027066LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA027066LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.