



ZNF98 Antibody, Biotin conjugated

Product Code	CSB-PA027106LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NK75
Immunogen	Recombinant Human Zinc finger protein 98 protein (1-135AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 98 (Zinc finger protein 739) (Zinc finger protein F7175), ZNF98, ZNF739
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF98
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.