



ompA Antibody, HRP conjugated

Product Code	CSB-PA028238LB01ZBR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	T2DH55
Immunogen	Recombinant Chlamydia pecorum Major outer membrane porin protein (24-374AA)
Raised In	Rabbit
Species Reactivity	Chlamydia pecorum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane porin, ompA
Immunogen Species	Chlamydia pecorum
Research Area	Others
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.