



# ompA Antibody, Biotin conjugated

Product Code	CSB-PA028238LD01ZBR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	T2DH55
Immunogen	Recombinant Chlamydia pecorum Major outer membrane porin protein (24-374AA)
Raised In	Rabbit
Species Reactivity	Chlamydia pecorum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane porin, ompA
Immunogen Species	Chlamydia pecorum
Research Area	Others
Target Names	ompA