



TMEM109 Antibody, Biotin conjugated

Product Code	CSB-PA03899D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BVC6
Immunogen	Recombinant Human Transmembrane protein 109 protein (34-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 109 (Mitsugumin-23) (Mg23), TMEM109
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMEM109
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.